



Model : AGIS Power Assist PS-HA-2



This power stroll assist can fit into most brand of Manual wheelchairs and it is easy to install and disconnect. It will take the strain away and make pushing the wheelchair easy on uphill ramps and uneven road. It increase the safety on down hill slope with its build in motor brake that allow the wheelchair to stop on slope and not worry about it rolling away.

Specifications

Power Stroll overall size (L x W x H)	950 X 670 X 840mm
Standard Wheelchair Model	KY 868
Seat Size	18 Inch
Wheelchair Weight	14KG
Power Stroll weight	19.2KG
Drive Wheel	200 X 50
Controller	PG S Drive 45A
Battery	12Ah×2
Charger	2A DC24V
Charging time , hours	> 4hrs from 50%
Motor	150W
MAX. speed	6.4kph
Cruising range	16KM
Weight capacity	135kg
Maximum climbing angle	< 6°

* All features and specifications are subject to change without notice