

MaxiCare 2-1 T07-P06A Tubular Alternating Cells Mattress(A-B alternating)
Pressure Relief Mattress



The MaxiCare 4.5 inch tubular overlay medical air mattress is an advanced multi-functional microcontroller programmed **7.0-14L Dual pump** dynamic mattress overlay system. The pump has 30 minutes **static mode and alternate mode with 4 selectable cycle time**. To ensure patient's safety, the pump has **audiovisual** low pressure, power failure, and service alarms.

Pump		Mattress	
Model	P06A	Model	4.5" Overlay T07
Dimension	345*115*185mm(L*W*H)	Dimension	L80" x W34" x H4.5"
Supply Voltage	230/50Hz	No of Cells	18
Alternating Type	2 tubes alternating	Cell Material	Nylon-TPU
Flowrate (litre)	7 - 14L	Top Cover Material	Nylon-PU
Control Panel	Dial + Digital Keypad	Mattress Base Material	4-way stretch Polyester-PU
Mode	Alternate & Static	Pillow Section	Yes - top 3 cells
Full Cycle Time (min)	Selectable 6/9/12/25 mins	Cell on Cell Design	No
Alarm	Audiovisual Alarms	Independent Low Air Lose	Yes
Press Ranges	10-56mmHg	CPR	Quick Release
Power Cord Type	Undetachable	Maximum Weight	145KG
Warranty	1 yr	Warranty	1 yr

The SuperCare Alternating enhances the comfort for the patient by fully re-inflate all sets of air cells before the next set of air cells deflates. Adjustable analog dial comfort control.

The innovative Low Air Loss design effectively reduces the moisture and skin heat on the support surface by periodically providing a high volume of low air loss relief.

Four Selectable 6, 9, 12 & 25 minutes alternating cycle time and static mode for nursing procedure.
Audiovisual alarms (power failure, low pressure, and service alarm).

Medical grade vapor permeable and spill resistance with fire retardant cover keeps the moisture away from patient's skin. Two high quality molded-tooth zippers cover for easy replacement and laundry.

Dual-way "Quick-CPR" provides immediate mattress deflation in case of emergency – FDA; CE; ISO13485